


Galt|STICK®



Embedded radiopaque marker band on all sheaths



Kits available with guidewire



Echogenic Trocar needle

The GaltStick® Introducer System is a 6 French Coaxial sheath assembly with stiffening stylet. The GaltStick® allows for the introduction of a .018" guidewire and .038" working guidewire in percutaneous procedures.

- Embedded radiopaque marker band on all sheaths for enhanced visibility
- Smooth dilator/sheath transition
- Variety of kit offerings
- Hydrophilic or silicone coated sheath configurations
- 6Fr outer sheath, 4Fr dilator, 20cm usable length
- Echogenic Trocar needle



For more information about the GaltStick® or to set up a meeting with a Galt Sales Consultant call (800) 639-2800.

Galt|STICK[®]



GaltStick[®] Introducer Set with Radiopaque Band

Model Number	Hydrophilic	.018" x 60cm SS/ Gold guidewire	.018" x 60cm NT/ Gold guidewire	.038" x 145 PTFE 3mm "J" Straight guidewire	Trocar Needle, echo	Quantity
KIT-030-04		x		x	x	5
KIT-030-05			x	x	x	5
KIT-030-06		x			x	10
KIT-030-03			x		x	10

Hydrophilic Coated GaltStick[®] Introducer Set with Radiopaque Band

Model Number	Hydrophilic	.018" x 60cm SS/ Gold guidewire	.018" x 60cm NT/ Gold guidewire	.038" x 145 PTFE 3mm "J" Straight guidewire	Trocar Needle, echo	Quantity
KIT-036-04	x	x		x	x	5
KIT-036-05	x		x	x	x	5
KIT-036-06	x	x			x	10
KIT-036-03	x		x		x	10

Accessories

Model Number	Description	Quantity
NDL-101-00	21ga x 15cm Trocar, Echogenic	5

